

LMCE KSLM

Laboratory
Medicine
Congress
& Exhibition

2016 & 57th ANNUAL MEETING

October 26-28, 2016 The K-Hotel, Seoul, Korea



R-Biopharm Korea

Product Highlights

Vision ESR analyzer



As a new generation of dynamic ESR analyzer, Vision can test up to 64 samples at the same time. With revolutionary breakthrough in technology, primary EDTA tubes can be directly used for determination and it only takes 20 minutes to obtain the results with have good correlation with the results of Westergren Method.

BIO SOLEA 2



BIO SOLEA 2 is semi-automated coagulation analyzer. BIO SOLIA 2 can be used with all reagent and work with any standard pipette. BIO SOLIA 2 is easy to operate(automatic incubation after insertion of samples cups, start upon reagent dispensing, calculated and displayed).